

# EyeMax 2

SENSING

ANALYSIS

DELIVERY

*"Optimum comfort and protection  
during phototherapy treatments."*





## EYEMAX 2 OVERVIEW

Our latex-free **EyeMax 2** provides optimum comfort and protection for patients during phototherapy treatments. The one-piece, wraparound design has two attachment points that can be adjusted independently; this prevents unwanted movement and allows for a perfect fit. The ergonomic eye-protection pad was specially designed to prevent ultraviolet light leakage.

## PART NUMBERS

EyeMax 2 Regular (20 pack) . . . . . R300P01  
EyeMax 2 Premie (20 pack) . . . . . R300P02  
EyeMax 2 Micro (20 pack) . . . . . R300P03



**1. Velcro Strapping**

*Adjustable band for optimum comfort*

**2. Eye Pads**

*Bigger, less pressure, and better UV protection*

**3. Latex Free**

*Helps to prevent irritant contact dermatitis*



**REGULAR** (12.6"-14.9")  
(Occipital Frontal Circumference)



**PREEMIE** (10.4"-12.6")  
(Occipital Frontal Circumference)



**MICRO** (7.87"-10.4")  
(Occipital Frontal Circumference)

**Specifications**

Regular	12.6" - 14.9" (Occipital Frontal Circumference)
Preemie	10.4" - 12.6" (Occipital Frontal Circumference)
Micro	7.87" - 10.4" (Occipital Frontal Circumference)