

MEDICAL AUTOMATIC ELECTRONIC BLOOD PRESSURE MONITOR

YE991



Large Backlit display



Body induction



Printing function

- Animated voice prompt, for automatic measurement
- Voice broadcast function, broadcasting the measured blood pressure in a timely manner
- Printer function, printing the measured blood pressure value
- Supporting multiple transmission modes
Intelligent pressurization and fixed pressurization dual-mode measurement

| | | | |
|----------------------------------|---------------------------------------------------------------|---------------------------------------|---|
| Applicable arm circumference | 17~45 cm | Voice broadcast function | ● |
| Pressure measuring range | 0~300mmHg(0~40kPa) | Printing function | ● |
| Systolic | 40~280mmHg | Blood pressure level display function | ● |
| Diastolic | 20~210mmHg | Automatic zeroing function | ● |
| Pulse measuring range | 40~240 time/min | Arm cylinder rotation function | ● |
| Accuracy of pressure measurement | ±2mmHg(±0.267kPa) | Cuff change prompt function | ● |
| Accuracy of pulse measurement | ±2%or±2 time/min | Elbow position prompt function | ● |
| Measurement time | About 35s | Average measuring function | ● |
| Data transmission | Bluetooth, Wi-Fi, USB drive, Ethernet, and RS-232 serial port | Proximity prompt function | ● |
| Display mode | LCD | Emergency stop function | ● |
| Systolic pressure display | ● | Indication for irregular heartbeats | ● |
| Diastolic pressure display | ● | Memory retrieval | ● |
| Pulse rate display | ● | | |

The parameters of YE991 are the same as the YE990, only there is no 'Prompt for movement error' function.